


PLENARIS

50 BODY HA

Crosslinked Hyaluronic Acid gel _____ 24mg/mL


Enhancing Body Volume and Contouring with Precision



PRODUCT INFORMATION

- **Package:** 1 Box (50 mL x 1 vial)
- **Hyaluronic Acid Concentration:** 24 mg/mL
- **Viscoelasticity:** 400 Pa.s – 700 Pa.s
- **Duration:** Over 12 months
- **Needle Compatibility:** 16G – 18G
- **Shelf-Life:** 24 months

High-Performance Viscoelasticity



With a viscosity range of 400 Pa.s – 700 Pa.s, PLENARIS 50 BODY HA ensures the filler remains precisely in the injected area. This targeted action minimizes the risk of product migration, delivering a natural and uniform result.

Premium-Grade Hyaluronic Acid

Formulated with pharmaceutical-grade, high-purity stabilized hyaluronic acid, originally designed for ophthalmology applications, PLENARIS 50 BODY HA ensures the utmost safety, compatibility, and reliability.

Extended Duration with Advanced Crosslinking

The high crosslinking rate prevents premature degradation by external factors, offering effects that last over 12 months. Patients enjoy enduring results with just one treatment session.



PLENARIS

50 BODY HA

Crosslinked Hyaluronic Acid gel _____ 24mg/mL

Enhancing Body Volume and Contouring with Precision

INDICATIONS

- **Volume Enhancement:** Breasts, buttocks, hips, etc.
- **Symmetry Restoration:** Asymmetric breasts, buttocks, calves
- **Liposuction Aftercare:** Subcutaneous tissue and volume restoration

BENEFITS

- **Non-surgical Volume Enhancement:** Perfect for patients seeking alternatives to invasive procedures for areas like breasts and buttocks.
- **Symmetry and Shape Correction:** Improves the aesthetic appearance of calves, hips, breasts, and other body areas.
- **Skin Rejuvenation:** Smooths skin, combats cellulite, and addresses unevenness after liposuction.
- **Quick Recovery:** Faster healing and downtime compared to traditional plastic surgery.

TREATMENT AREAS

